

# Disposable Instruments

for use with PrimeSight™ Endoscopy Systems




## Biopsy Forceps

Cogentix Medical Single-use Biopsy Forcep is intended for use in collecting tissue samples in histological exams.

- » Surgically sharp edge
- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9007	1.7mm	160cm	Endoscopy	Box of 10




## Tripod Grasper

The Cogentix Medical Retrieval Forcep is intended for use in removing foreign debris.

- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9008	1.7mm	150cm	Endoscopy	Box of 10



## Retrieval Baskets

Cogentix Medical Retrieval Baskets are intended for use in removing foreign debris.

- » Hexagonal and elliptical designs available
- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9009 (Elliptical)	1.7mm	150cm	Endoscopy	Box of 10
07-9010 (Hexagonal)	1.7mm	150cm	Endoscopy	Box of 10

For more information about PrimeSight Endoscopy Systems  
visit [primesightendoscopy.com](http://primesightendoscopy.com)





## Aspirating Syringe

The Cogentix Medical Locking Aspirating Syringe is designed with a unique locking system for maintaining negative pressure.

- » Universal Luer-lock fitting
- » Disposable to prevent cross contamination
- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9003	N/A	N/A	Endoscopy	Box of 5



## Cytology Brushes

Cogentix Medical Cytology Brushes offer the endoscopist the ideal combination of high yield and ease of use. The unique design of the brush allows maximum specimen yield with minimum cellular deposits left in the sheath.

- » Teflon sheathing for ease of insertion and manipulation during procedures
- » Tough stainless steel, high tempered wire to prevent kinking
- » Radiopaque ball tip is designed to reduce mucosal trauma
- » Thumb ring for easy operations
- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9001	1.7mm	80cm	Gastroenterology Pulmonology	Box of 10
07-9006	1.7mm	10cm	ENT (Nasal)	Box of 10



## Injection Needle

The Cogentix Medical Injection Needle features a uniquely designed extension and retraction system, ensuring accurate and quick release at the injection site.

- » Accurate, dependable full needle extension
- » Single hand operation
- » Protective radiopaque tip eliminates risk of scope damage
- » Packaged sterile for single use

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9005	1.8mm	80cm	Gastroenterology	Single



## Transbronchial Aspiration Needle

The Cogentix Medical Transbronchial Aspirating Needle is designed for maximum specimen retrieval when used in conjunction with our self-locking Aspirating Syringe.

- » 21 gauge needle for maximum specimen retrieval
- » High specimen yield when negative pressure is maintained with Aspirating Syringe
- » Stop feature for accurate extension and retraction

Catalog Number	Working Diameter	Working Length	Specialty	Quantity
07-9002	1.8mm	80cm	Pulmonology	Single

Connect with Cogentix Medical



Rx Only. PrimeSight is a trademark of Cogentix Medical. © 2015 Cogentix Medical. All rights reserved.  
70021A 10/15

**Cogentix**  
Medical

Cogentix Medical  
5420 Feltl Road  
Minnetonka, MN 55343

cogentixmedical.com

TEL 866.258.2182  
FAX 866.255.4522