




Sterile Endoscopy Solution for Every Patient



The sterile, single-use EndoSheath protective barrier includes the channel, ports and seals and is used with the D-shaped PrimeSight™ Endoscope.



Prior to each procedure the PrimeSight Endoscope is inserted into the EndoSheath.



The EndoSheath completely isolates the PrimeSight Endoscope from the patient during procedure.

Proven Performance. Simplified Safety.

- » FDA-cleared microbial barrier for flexible endoscopy
- » Eliminates the need for HLD between procedures, reducing the time and complexity of reprocessing
- » Over 5 million used without a single reported complaint of patient-to-patient cross contamination
- » Proven effective barrier against organisms as small as 27 nanometers per FDA guidance and testing
- » 100% leak tested during manufacturing with a 0% acceptable failure rate
- » 30 plus research studies performed on EndoSheath since 1992

Recognition in National Standards and Guidelines

- » ANSI/AAMI ST91:2015 National Standards for Flexible and Semi-Rigid Endoscope Processing in Health Care Facilities
- » HICPAC Guideline for Disinfection and Sterilization in Healthcare Facilities
- » AUA/SUNA Guidelines for Cystoscopy Reprocessing

The EndoSheath Protective Barrier is available for PrimeSight Systems in Bronchoscopy, Cystoscopy, Esophagoscopy, and Laryngoscopy

PrimeSightEndoscopy.com

Cogentix
Medical

Cogentix Medical
5420 Feltl Road
Minnetonka, MN 55343

cogentixmedical.com

TEL 866.258.2182
FAX 866.255.4522