




## Sterile Endoscopy Solution for Every Patient



The sterile, single-use EndoSheath protective barrier includes the channel, ports and seals and is used with the D-shaped PrimeSight™ Endoscope.



Prior to each procedure the PrimeSight Endoscope is inserted into the EndoSheath.



The EndoSheath completely isolates the PrimeSight Endoscope from the patient during procedure.

### Proven Performance. Simplified Safety.

- » FDA-cleared microbial barrier for flexible endoscopy
- » Eliminates the need for HLD between procedures, reducing the time and complexity of reprocessing
- » Over 5 million used without a single reported complaint of patient-to-patient cross contamination
- » Proven effective barrier against organisms as small as 27 nanometers per FDA guidance and testing
- » 100% leak tested during manufacturing with a 0% acceptable failure rate
- » 30 plus research studies performed on EndoSheath since 1992

### Recognition in National Standards and Guidelines

- » ANSI/AAMI ST91:2015 National Standards for Flexible and Semi-Rigid Endoscope Processing in Health Care Facilities
- » HICPAC Guideline for Disinfection and Sterilization in Healthcare Facilities
- » AUA/SUNA Guidelines for Cystoscopy Reprocessing

The EndoSheath Protective Barrier is available for PrimeSight Systems in Bronchoscopy, Cystoscopy, Esophagoscopy, and Laryngoscopy

[PrimeSightEndoscopy.com](http://PrimeSightEndoscopy.com)

**Cogentix**  
Medical

Cogentix Medical  
5420 Feltl Road  
Minnetonka, MN 55343

[cogentixmedical.com](http://cogentixmedical.com)

TEL 866.258.2182  
FAX 866.255.4522