


COGENTIX MEDICAL

# *Stone Extractors*


## Reliable Stone Extraction

Cogentix Medical offers a variety of stone extractors to aid in stone management for use in ureteroscopy procedures.




**Cogentix**  
Medical

## Stone Extractors



### Tipless Nitinol Baskets

Available in 1.3 Fr, 1.9 Fr, and 2.2 Fr, the tipless baskets are ideal for the removal of upper urinary tract and renal stones. The small size of the baskets improve irrigation and visibility. The tipless configuration allows for close stone proximity for less traumatic retrieval of calyceal calculi. The flexible shaft allows for enhanced scope deflection.



### Flat and Helical Nitinol Baskets

Available in 2.5 Fr, these baskets are designed for rigid and semi-rigid ureteroscopy procedures. Both the flat wire and the helical design feature ergonomic, removable handles. Choose the wire design that best meets your stone manipulation needs.

## Product Information

CATALOG NUMBER	DESCRIPTION
Tipless Nitinol Stone Extractors	
09-0023	Nitinol stone extractor, 12 mm tipless basket, 4 wires, 1.3 Fr
09-0024	Nitinol stone extractor, 16 mm tipless basket, 4 wires, 1.9 Fr
09-0025	Nitinol stone extractor, 12 mm tipless basket, 4 wires, 1.9 Fr
09-0026	Nitinol stone extractor, 12 mm tipless basket, 4 wires, 1.9 Fr, criterium handle
09-0028	Nitinol stone extractor, 16 mm tipless basket, 4 wires, 2.2 Fr
Flat and Helical Nitinol Stone Extractors	
09-0019	Nitinol stone extractor, 16 mm straight basket, 4 wires, 2.5 Fr
09-0021	Nitinol stone extractor, 16 mm helical basket, 4 wires, 2.5 Fr